## **Cancer Screening Tip of the Month – May 2025**

Toronto Central Regional Cancer Program, Cancer Screening eBulletin

## Alert: HPV Collection Vials in the Ontario Cervical Screening Program (OCSP)

- ThinPrep® Pap Test vials (white and purple vials) should be used for all testing in the OCSP, including provider-collected cervical specimens and self-collected vaginal specimens at OCSP self-collection sites.
- ThinPrep® PreservCyt vials (green vials) are not approved by Health Canada for the transportation and preservation of gynecologic samples (e.g., cervical and vaginal samples).

ThinPrep® Pap Test collection vials

✓ Self-collected vaginal specimens
✓ Provider-collected cervical specimens

ThinPrep® PreservCyt collection vials

X Should NOT be used for any testing within the OCSP

Click here to access HPV testing tools and resources for primary care providers.